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What is a mixSorb doing?

The mixSorb is a dynamic sorption analyzer to investigate adsorbent materials for gas mixture or gas/vapor mixture separation processes.

- investigations to technical relevant parameters (sorption capacity, selectivity, sorption kinetics, regenerability, cycle stability ...)
- technical processes in “Labscale”

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release:
Dec, 2015

mixSorb L & LHP

**large amount of sample/granulate (~ 130 ml)
application-near measurements**

→ companies/producers and institutes for applied sciences

release:
Oct, 2018


mixSorb S & mixSorb SHP

**small amount of sample/powder (~ 1 ml)
basic research**


→ institutes for basic research and R&D for powders

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
Different segments of a breakthrough curve

- Determination of technical usable sorption capacity
- Can be used as benchmark for separation performance of adsorbents


- Mass Transfer coefficient, axial dispersion, shape of isotherm
- Heat effects, heat dissipation
- The time interval of mass transfer zone has to be minimized

- Determination of saturation capacity
- By assuming of thermodynamic controlled system → Measurement of isotherms possible

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
simple breakthrough curve

- binary mixture:
 - CO₂/He (non-adsorbing carrier gas)
 - Pure component equilibria+ kinetic information
 - CO₂/CH₄ (adsorbing carrier gas)
 - Preloading of sample with pure CH₄
 - Incomplete Determination of the system (evaluation mostly simple)
 - Partial loading for CO₂ (mixture sorption data)

breakthrough curves with displacement


- ternary mixture:
 - CO₂/CH₄/He (non-adsorbing carrier gas)
 - Displacement of less adsorbed component
 - Partial Desorption, Role-Up Effects
 - Complete Determination (evaluation complex)
 - CO₂/CH₄/N₂ (adsorbing carrier gas)
 - Preloading of sample with pure N₂
 - Partial loading for CO₂ and CH₄
 - Incomplete ternary mixture data

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
Different experiments with mixSorb L and their results:


Breakthrough time, Mass Transfer


Competitive Adsorption, Displacement

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




- Flow through the sample with a predefined gas mixture
- Measurement of equilibrium data at a specified pressure and gas composition

➔ **Advantage over static volumetry, no precalculations necessary**

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
Table 1 Summary of all conducted calibration experiments

#	Adsorber Filling	Gas	T / °C	Flow rate / mL min ⁻¹	Purpose
1	empty	He	25	20	▪ Determination of calibration factor He ▪ Manifold volume
2	empty	He	50	20	
3	Cu beads	He	25	20	
4	empty	CH ₄	25	3	Conversion factor over the whole pressure range
5	empty	CO ₂	25	3	Conversion factor over the whole pressure range

Dynamic Manometric Physisorption: Principle

- **Continuous** pressure increase
- Possible with mixSorb S and SHP by default
- Flow slow enough to equilibrate → 5 mbar/s appears to be okay

$$\frac{\dot{V}}{V_{\text{free}}} \cdot 1000 \text{ mbar}$$

Calibrations

- \dot{V} : exact flow rates for Helium and the adsorptive
→ needs to be performed for each experimental setup
- $V_{\text{warm fs}}$ for Manifold
→ needs to determine once per instrument

Corrections


- Free space determination with Helium to correct free cold space $V_{\text{cold fs}}$

Considerations

- Apply good **EOS**

Dynamic Manometric Physisorption: Pressure Rise


Example: CO₂ Isotherm

Step 1: Data recording

- For every datapoint:
 - Time
 - Pressure
 - Flow (@ STP)
 - Manifold temperature („warm“)
 - Sample temperature („cold“)
- Apply EOS to calculate densities
- $\rho_{p,warm} \cdot V_{free,warm} + \rho_{p,cold} \cdot V_{free,cold} = n_{detected}$
- $\int_0^t \dot{V}_t dt \cdot \rho_{STP} = n_{dosed}$

Dynamic Manometric Physisorption: Amount balance


Example: CO₂ Isotherm

Step 2: calculate amounts

- For every datapoint we can now add:
 - Detected amount
 - Dosed amount up to time t
- Draw detected amount over dosed amount
- When **no adsorption** takes place: straight line with slope of „1” → correction of non-idealities

Dynamic Manometric Physisorption: Isotherm

Example: CO₂ Isotherm


Step 3: Build isotherm

- $n_{\text{ads}} = n_{\text{dosed}} - n_{\text{detected}}$
- Divide by sample mass
- Reduction of cold free space by adsorbate can also be incorporated

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Dynamic Manometric Physisorption: Comparison

Comparison to gravimetric sorption isotherms:

- Promising
- Very good agreement between 1 and 10 bar
- Measurement time **5 h** vs. 20 h with gravimetry

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
Leipzig

Adsorption Week

Adsorption and Characterization of porous solids

Keynote Speaker:
Stefan Kaskel (University
Dresden, Germany)


- Standardizace: Adsorption information file
- Daniel W. Siderius (NIST) – NIST Data Resources for Adsorption
- Stefan Kaskel (TU Dresden)
 - Paralelní charakterizace adsorpce s XRD, NMR, XAS, EPR, UV/Vis/Raman
 - Negativní adsorpce (kontrakce struktury)
- Jean Rouquerol (Aix Marseille Université, Lab Madirel)
 - Constant rate thermal analysis – thermal treatment

Adsorption Information File - AIF (Stefan Kaskel)

<https://adsorptioninformationformat.com/>

Adsorption information can be shared using the STAR format, which is used by the crystallographic information file (CIF).

```
# comment goes here                                     comment
data_raw2aifv003                                         experimental
_exptl_operator 'Nicole'
_exptl_date 2009-09-17T19:56:00
_exptl_instrument 'Autosorb Station 1'
_exptl_adsorptive Nitrogen
_exptl_temperature 77.3
_exptl_sample_mass 0.0339
_sample_id 'nk_DUT-6_LP_N2_114pkt'                       sample
_sample_material_id 'DUT-6'
_units_temperature K                                     units
_units_pressure Pa
_units_mass g
_units_loading mmol/g
loop_                                                     adsorption
  _adsorp_pressure
  _adsorp_p0
  _adsorp_amount
0.269367243408 101860.98004799998 0.006484305926579284
0.27275763204 101860.98004799998 0.008217577878132952
0.273485572344 101860.98004799998 0.009928476539499517
loop_                                                     desorption
  _desorp_pressure
  _desorp_p0
  _desorp_amount
97113.50270639999 101860.98004799998 57.919842028778376
90734.559156 101860.98004799998 57.83008565285129
85586.58132479999 101860.98004799998 57.74993663094573
```

[Adsorption Information File - Metal-Organic Frameworks - International Commission \(mof-international.org\)](https://www.mof-international.org/)

Daniel W. Siderius (NIST): NIST Data Resources for Adsorption

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[HTTPS://ADSORPTION.NIST.GOV](https://adsorption.nist.gov)

Daniel W. Siderius (NIST): NIST Data Resources for Adsorption

Further Reading

Park, Howe, Sholl: **“How Reproducible Are Isotherm Measurements in Metal–Organic Frameworks?”**

Chem Mater, 2017 (10.1021/acs.chemmater.7b04287)

Bingel, Chen, Agrawal, Sholl: **“Experimentally Verified Alcohol Adsorption Isotherms in Nanoporous Materials from Literature Meta-Analysis”**

J Chem Eng Data, 2020 (10.1021/acs.jced.0c00598)

Iacomi and Llewellyn: **“pyGAPS: a Python-based framework for adsorption isotherm processing and material characterization”**

Adsorption, 2019 (10.1007/s10450-019-00168-5)

Evans et al. **“A Universal Standard Archive File for Adsorption Data”**

Langmuir 2021 (10.1021/acs.langmuir.1c00122)


Bucior et al. **“Identification schemes for metal–organic frameworks to enable rapid search and cheminformatics analysis”**

Cryst Growth Des, 2019 (10.1021/acs.cgd.9b01050)

[HTTPS://ADSORPTION.NIST.GOV](https://adsorption.nist.gov)

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111 K,
boiling point of CH_4

- Desorption step in Adsorption branch (163 cm³/g!)
- $dn/dp < 0$

S. Krause, V. Bon, F.-X. Coudert, G. Maurin, S. Kaskel et. al. *Nature* **2016**, 532, 348–352.

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V. Bon, E. Brunner, A. Pöpl, S. Kaskel,
“Unraveling Structure and Dynamics in Porous Frameworks via Advanced In Situ Characterization Techniques”.
Adv. Funct. Mater. **2020**, 1907847. <https://doi.org/10.1002/adfm.201907847>

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ROUQUEROL F., ROUQUEROL J., SING K.S.W., LLEWELLYN P., MAURIN G. « Adsorption by powders and porous solids: principles, methodology and applications” 2nd edition , Academic Press, 2014, 626 p.

SORENSEN T. O., ROUQUEROL J. Eds. « Sample Controlled Thermal Analysis (SCTA): principle, origins, goals, multiple forms, applications and future” Kluwer Academic Publishers, Dordrecht, 2003, 252 p

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Gas adsorption measurements: All the things that can go wrong...

<https://www.3p-instruments.com/free-webinars-gas-adsorption-measurements/>

Ing. Michal Dudák, Ph.D.

ANAMET s.r.o.

dudak@anamet.cz